



22023.03.22 Ver.3

Patient Transfer System



Instruction for use

EZ-900B/910B Patient Handling Belt















Important:

Read these instructions carefully before using the EZ-900B/910B Patient Handling Belt.

Manufacturer

EZ-Care Solution, Co., Ltd.
12F, No. 171, Sec. 4, Nanjing E. Rd.,
Songshan Dist., Taipei City 105, Taiwan
TEL: +886-2-2659-2277; FAX:+886-2-2798-8321;
E-mail: sales@ezcaresolution.com
www.ezcaresolution.com

All of the aids produced by EZ-Care Solution Co., Ltd fulfill the following norms for the medical products Class 1 (MDR (EU) 2017/145).

Explanation of symbols in Instruction for use and on labels			
	Model Number		Warning
	Lot No.		Hang Dry
	Read the user manual		Manufacturer
	No sharps		European Representative
	Do not bleach		Medical device
	Do not iron		CE-marking
	Machine washable		Unique Device Identification

Product Label

EZ-Go®
Patient Transfer System

Patient Handling Belt

REF EZ-900B

LOT 230301

 Do Not high temp sterilize

Textile: Nylon
Inner Board: EPE

MADE IN TAIWAN

 **EZ-CARE**

EZ-Care Solution Co., Ltd.
12F, No.171, Sec 4, Nanjing
E. Rd., Songshang Dist.,
Taipei Taiwan
www.ezcaresolution.com
sales@ezcaresolution.com

EC REP Obelis s. a.
Boulevard Général Wahis,
53 1030 Brussels, Belgium
Tel: +32 (2) 73 25 954
sales@obelis.net

MD **CE**

UDI


471704851EZ-9004C

EZ-Go®
Patient Transfer System

Patient Handling Belt

REF EZ-910B

LOT 230301

 Do Not high temp sterilize

Textile: Nylon
Inner Board: EPE

MADE IN TAIWAN

 **EZ-CARE**

EZ-Care Solution Co., Ltd.
12F, No.171, Sec 4, Nanjing
E. Rd., Songshang Dist.,
Taipei Taiwan
www.ezcaresolution.com
sales@ezcaresolution.com

EC REP Obelis s. a.
Boulevard Général Wahis,
53 1030 Brussels, Belgium
Tel: +32 (2) 73 25 954
sales@obelis.net

MD **CE**

UDI


471704851EZ-9004C

Contents

1. Introduction	5
1.1 Intended Use	5
1.2 Intended Environment	5
1.3 Safety	5
1.4 Warranty	5
1.5 Product Structure	6
1.6 Product Specification	6
2. Operation	7
2.1 Placement	7
2.2 Transfer from Bed to Chair.....	8
2.3 Sitting to Standing.....	9
2.4 Walking Training	9
3. Cleaning Instruction	10
4. Maintenance	10

1. Introduction

1.1 Intended Use

EZ-900B/910B is a durable patient transfer handling belt that helps caregivers assist patients (1) from sitting to standing, (2) transferring from bed to chair, and (3) walking training. The attached handles enable caregivers to hold the patient without the need to grasp their clothes. EZ-900B/910B provides safer and easier patient handling and helps reduce work-related injuries for caregivers.

1.2 Intended Environment

Acute Care, Long-term Care, Home Care

1.3 Safety



A failure to comply with the following safety points and further definitions in these instructions for use could lead to dangerous situations. Please read carefully.



Please consult with the patient's care team in advance to ensure the product can be used, and in combination with other devices, that is safe for the patient and the caregivers.



Always complete a patient risk assessment prior to all patient-handling activities.



It's important that caregivers receive manual handling training in advance.



The product is not designed for lifting patients.



Ensure the handling belt is securely fastened before transferring.




When transferring to another device (such as from bed to wheelchair), ensure there is no gap between two surfaces and two surfaces are leveled.

1.4 Warranty

All EZ-GO products are supplied with a 1-year warranty, on construction and/or material faults unless otherwise stated. In the event of varying warranty terms, the warranty starts from the invoice date. The warranty does not cover failures and damage not caused by the product's defects such as improper cleaning.

1.5 Product Structure



1.6 Product Specifications

Model No./ Specification ($\pm 3\%$)	EZ-900B	EZ-910B
Size	107 x 16 cm	128 x 16 cm
Net Weight	610 g	680 g
Suitable Waist Size	26 ~35 inches	35 ~43 inches

Material:

Textile: Nylon Anti-slip mesh: Eco PER


Buckle: POM Inner Board: EPE board

Lifetime: under normal storage and use condition, 2 years.


2. Operation

2.1 Placement

- Place EZ-900B/910B around the patient's waist. Adjust the belt length to fit the patient and securely fasten the belt. Attach one end of the crotch strap to the main belt at the back. Have the crotch strap go through the inner thighs and attach the other end of the crotch strap to the front of the main belt. Adjust the strap length to fit the patient and securely fasten the strap. Apply the same steps to the other strap.




2.2 Transfer from bed to chair

 The patient must be cooperative and predictable, able to bear weight on both legs and take small steps.

 Avoid lifting patients.

- Lower the bed and ensure that brakes are applied. Also, ensure brakes are applied on the wheelchair.
- Place the wheelchair next to the bed at a 45-degree angle and apply brakes. If the patient has weakness on one side, place the wheelchair on the strong side.
- Sit the patient on the side of the bed with his or her feet on the floor. The patient's feet should be in between the caregiver's feet.
- Place hands on the waist to assist into a standing position.
- As the patient leans forward, grasp the handling belt on the side of the patient, with the arms outside the patient's arms. Position the legs on the outside of the patient's legs.
- Using a rocking motion, help the patient stand by shifting weight from the front foot to the back foot, keeping elbows in and back straight.
- Once standing, have the patient take a few steps back until they can feel the wheelchair on the back of their legs. Have the patient grasp the armrest of the wheelchair and sit down.



2.3 Sitting to Standing

- Assist the patient to move forward in a chair/bed, using an appropriate method.
- Assist the patient into a “nose over toes” position, with both feet firmly on the floor, one slightly in front of the other.
- Caregivers firm the feet on the floor with the feet slightly apart and the outer foot ahead. Flex the knees and hips to lower the position. Grasp the appropriate handle on the belt, with one hand. DO NOT put the whole hand through the loop.
- Give clear instructions to the patient. A gentle rocking motion can be used.
- Both caregivers move together in the direction of the movement, first forward and then upwards.



Ensure the patient is able to provide the majority of standing effort before use.



Ensure the caregivers have had proper safe patient handling training.




Stabilize the move before proceeding into a walk. If the user is unsteady, have the patient sit down.

2.4 Walking training

- Caregivers firm the feet on the floor with the feet slightly apart and the outer foot ahead, ready to step with the patient.
- Grasp the appropriate handle on the belt, with one hand. DO NOT put your whole hand through the loop.
- To initiate walking, keep as close to the patient as necessary, usually slightly behind, with the hip close to or touching the patient’s pelvis, as you walk forward. It is advisable to give clear instructions to the patient throughout the walking operation.



If the patient begins to fall from a standing position, do not attempt to stop the fall or catch the patient. Instead, control the fall by lowering the patient to the floor.



3. Cleaning Instruction

- Machine wash with a temperature up to 30°C.
- Do not bleach.
- Do not iron.
- Do not high-temperature sterilization.



4. Maintenance

- Store at room temperature.
- EZ-900/910 must be checked periodically.
 - ✓ Visual Check: Ensure the materials, seams, and stitching are intact/faultless.
 - ✓ Test Handles: Draw hard in the handles in opposite directions to ensure all handles are well attached.